

Product Datasheet

Phencyclidine (PCP) Antibody P01-99-63M-P

Product Name	Phencyclidine (PCP) Antibody
Catalog Number	P01-99-63M-P
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Phencyclidine-BSA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG3
Clonality	Monoclonal
Buffer	0.05 M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Specificity	Phencyclidine
Matched Pairs	P81-99-61A (capture) or P91-99-61H (detection)
Preservatives	0.05% Sodium Azide