

Product Datasheet

Phenobarbital Antibody P01-99-81M-P

Product Name	Phenobarbital Antibody
Catalog Number	P01-99-81M-P
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Phenobarbital-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Buffer	0.05 M Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	Phenobarbital
Matched Pairs	P81-99-81A (capture) or P91-99-81H (detection)
Preservatives	0.05% Sodium Azide