

Product Datasheet

Anti-Human CD126 Monoclonal Antibody, Azide Free Clone B-R6 852.030.005

Product Name	Anti-Human CD126 Monoclonal Antibody, Azide Free Clone B-R6
Catalog Number	852.030.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, IP, Functional assay, IHC, FC
Immunogen	gp80 transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Interleukin-6 Ra subunit, IL-6 Ra, IL-6 R, a 80 kDa protein
Biological Activity	Inhibits IL-6 induced proliferation of XG-1 cells
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells