

Product Datasheet

Carisoprodol (Soma) Antibody P02-99-73R-IF

Product Name	Carisoprodol (Soma) Antibody
Catalog Number	P02-99-73R-IF
Description	Protein A Purified Polyclonal Antibody, Ig Fraction
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Carisoprodol - BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05 M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	meprobamate 12%
Specificity	Carisoprodol
Preservatives	0.05% Sodium Azide