

Product Datasheet

Anti-Human CD126 Monoclonal Antibody, FITC Conjugated Clone B-R6 852.031.010

Product Name	Anti-Human CD126 Monoclonal Antibody, FITC Conjugated Clone B-R6
Catalog Number	852.031.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	gp80 transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5 % BSA and 0.09 % sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Interleukin-6 Ra subunit, IL-6 Ra, IL-6 R, a 80 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells