

## Product Datasheet

### Methaqualone Antibody P02-99-83R-IF

<b>Product Name</b>	Methaqualone Antibody
<b>Catalog Number</b>	P02-99-83R-IF
<b>Description</b>	Protein A Purified Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Methaqualone-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Polyclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
<b>Specificity</b>	Methaqualone
<b>Matched Pairs</b>	P82-99-81A (capture)
<b>Preservatives</b>	0.05% Sodium Azide