

Product Datasheet

Methylphenidate Antibody P03-99-03R-IF

Product Name	Methylphenidate Antibody
Catalog Number	P03-99-03R-IF
Description	Protein A Purified Polyclonal Antibody
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Methylphenidate-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Specificity	Methylphenidate
Matched Pairs	P83-99-01A (capture)
Preservatives	0.05% Sodium Azide