

Product Datasheet

Herpes Simplex Virus (HSV-gB) Antibody P1105

Product Name	Herpes Simplex Virus (HSV-gB) Antibody
Catalog Number	P1105
Description	Protein G Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤ -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gB