

Product Datasheet

Herpes Simplex Virus 2 (HSV 2-gC) Antibody P1111

Product Name	Herpes Simplex Virus 2 (HSV 2-gC) Antibody
Catalog Number	P1111
Species/Host	Mouse
Tested Applications	IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gC-2