

Product Datasheet

Herpes Simplex Virus 2 (HSV 2-ICP27) Antibody P1113

Product Name Herpes Simplex Virus 2 (HSV 2-ICP27) Antibody

Catalog Number P1113

Description Protein G Purified Monoclonal Antibody

Species/Host Mouse

Tested Applications WB, IFA

Brand EastCoast Bio

Form/Appearance Protein G Purified Monoclonal Antibody

Storage ≤ -20 °C. Avoid repeated freezing and thawing.

Isotype IgG

Clonality Monoclonal

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Specificity ICP27