

Product Datasheet

Cytomegalovirus (CMV-ICP36) Antibody P1202-1

Product Name	Cytomegalovirus (CMV-ICP36) Antibody
Catalog Number	P1202-1
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	$\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	ICP36