

# Product Datasheet

## Cytomegalovirus (CMV-ICP22) Antibody P1210

|                     |                                              |
|---------------------|----------------------------------------------|
| Product Name        | Cytomegalovirus (CMV-ICP22) Antibody         |
| Catalog Number      | P1210                                        |
| Species/Host        | Mouse                                        |
| Tested Applications | WB, IFA                                      |
| Brand               | EastCoast Bio                                |
| Form/Appearance     | Protein G Purified Monoclonal Antibody       |
| Storage             | ≤-20°C. Avoid repeated freezing and thawing. |
| Isotype             | IgG                                          |
| Clonality           | Monoclonal                                   |
| Buffer              | 10 mM Phosphate Buffered Saline, pH 7.4      |
| Specificity         | ICP22                                        |