

Product Datasheet

Anti-Human CD130 Monoclonal Antibody, FITC Conjugated Clone B-R3 852.061.010

Product Name	Anti-Human CD130 Monoclonal Antibody, FITC Conjugated Clone B-R3
Catalog Number	852.061.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	Natural soluble gp130
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5 % BSA and 0.09 % sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells