

## Product Datasheet

### Anti-Human CD130 Monoclonal Antibody, Unconjugated Clone B-R3 852.063.020

<b>Product Name</b>	Anti-Human CD130 Monoclonal Antibody, Unconjugated Clone B-R3
<b>Catalog Number</b>	852.063.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC, IHC
<b>Immunogen</b>	Natural soluble gp130
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Note</b>	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
<b>Isotype</b>	Mouse IgG2a Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognises Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein
<b>Hybridoma</b>	Myeloma X63/AG.8653 x Balb/c spleen cells