
Product Datasheet

Cotinine-4-BSA Antigen P82-99-02A

Product Name	Cotinine-4-BSA Antigen
Catalog Number	P82-99-02A
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	Distilled Water
Matched Pairs	PP194IF, PP195IF or PP198IF (detection)
Preservatives	None