

Product Datasheet

Methaqualone-BSA Antigen P82-99-81A

Product Name	Methaqualone-BSA Antigen
Catalog Number	P82-99-81A
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	Distilled Water
Matched Pairs	P02-99-81R-IF, 82R-IF or 83R-IF (detection)
Preservatives	None