

Product Datasheet

Cannabinoid (THC) Antibody PP123IF

| | |
|---------------------|---|
| Product Name | Cannabinoid (THC) Antibody |
| Catalog Number | PP123IF |
| Species/Host | Sheep |
| Tested Applications | LF, ELISA |
| Brand | EastCoast Bio |
| Form/Appearance | DEAE Purified Polyclonal Antibody; Ig Fraction |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Clonality | Polyclonal |
| Buffer | 20 mM Phosphate Buffer, pH 7.1 with 150 mM Sodium Chloride |
| Specificity | THC |
| Matched Pairs | LA020 (capture) |
| Preservatives | 0.05% Sodium Azide |