

# Product Datasheet

## Cannabinoid (THC) Antibody PP123IF

Product Name	Cannabinoid (THC) Antibody
Catalog Number	PP123IF
Species/Host	Sheep
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Polyclonal Antibody; Ig Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	20 mM Phosphate Buffer, pH 7.1 with 150 mM Sodium Chloride
Specificity	THC
Matched Pairs	LA020 (capture)
Preservatives	0.05% Sodium Azide