

# Product Datasheet

## Cannabinoid (THC) Antibody PP123IF

|                     |   |
|---------------------|---|
| Product Name        | Cannabinoid (THC) Antibody  |
| Catalog Number      | PP123IF   |
| Species/Host        | Sheep   |
| Tested Applications | LF, ELISA   |
| Brand               | EastCoast Bio   |
| Form/Appearance     | DEAE Purified Polyclonal Antibody; Ig Fraction  |
| Storage             | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Clonality           | Polyclonal  |
| Buffer              | 20 mM Phosphate Buffer, pH 7.1 with 150 mM Sodium Chloride                            |
| Specificity         | THC   |
| Matched Pairs       | LA020 (capture)   |
| Preservatives       | 0.05% Sodium Azide  |