

Product Datasheet

Cannabinoid (THC) Antibody PP123IF

Product Name	Cannabinoid (THC) Antibody
Catalog Number	PP123IF
Description	DEAE Purified Polyclonal Antibody; Ig Fraction
Species/Host	Sheep
Tested Applications	ELISA, LF
Alternative Names	Tetrahydrocannabinol
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Polyclonal Antibody; Ig Fraction
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	20 mM Phosphate Buffer, pH 7.1 with 150 mM Sodium Chloride
Specificity	THC
Matched Pairs	LA020 (capture)
Preservatives	0.05% Sodium Azide