

Product Datasheet

Anti-Human CD130 Monoclonal Antibody, Azide Free Clone B-K11 852.100.005

Product Name Anti-Human CD130 Monoclonal Antibody, Azide Free Clone B-K11

Catalog Number 852.100.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, IP, FC

Immunogen Natural soluble gp130

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognises the Gp130, common subunit for IL-6, IL-11, OSM, LIF, CNTF, CT-1 receptors, a 130-140 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells