

Product Datasheet

Cannabinoid (THC) Antibody PP124IF

Product Name	Cannabinoid (THC) Antibody
Catalog Number	PP124IF
Description	Purified Polyclonal Antibody
Species/Host	Sheep
Tested Applications	ELISA, LF
Immunogen	? 8 THC-BgG
Alternative Names	Tetrahydrocannabinol
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	50 mM Phosphate, pH 7.15
Matched Pairs	LA020 (capture)
Preservatives	0.02% Sodium Azide