

Product Datasheet

Cotinine-4 Antibody PP194IF

Product Name	Cotinine-4 Antibody
Catalog Number	PP194IF
Description	Protein A Purified Polyclonal Antibody, Ig Fraction
Species/Host	Rabbit
Tested Applications	ELISA, LF
Immunogen	Cotinine-4-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	P82-99-02A (capture) or P92-99-02H (detection)
Preservatives	0.05% Sodium Azide