
Product Datasheet

Ferritin, Polyclonal Antibody Anti-Human Affinity Purified RF-80A-10

Product Name	Ferritin, Polyclonal Antibody Anti-Human Affinity Purified
Catalog Number	RF-80A-10
Species/Host	Rabbit
Immunogen	High purity native human liver Ferritin
Brand	Lee Biosolutions
Form/Appearance	Liquid
Storage	2-8° C
Isotype	IgG
Purity	Affinity Purified
Source	Polyclonal Antibody
Shelf Life	2 years