

Product Datasheet

Anti-Human CD213a2 Monoclonal Antibody,Azide Free Clone B-D13 852.120.005

Product Name	Anti-Human CD213a2 Monoclonal Antibody,Azide Free Clone B-D13
Catalog Number	852.120.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, IHC, FC
Immunogen	IL-13Ra2 transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the IL-13 receptor a2 (IL-13 Ra2), a 44 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells