

Product Datasheet

dNTP Sets (100 mM solution) MX4501

Product Name	dNTP Sets (100 mM solution)
Catalog Number	MX4501
Description	MedixMDx deoxynucleotide triphosphate (dNTP) Set contains separately packaged 2'-deoxyadenosine-5'- triphosphate (dATP), 2'-deoxycytidine-5'-triphosphate (dCTP), 2'-deoxyguanosine-5'-triphosphate (dGTP), and 2'-deoxythymidine-5'-triphosphate (dTTP). The dNTP Set is a clear aqueous solution (pH 8.5) with greater than 99% purity. The solution is free of human and bacterial DNA as well as RNases and DNases. This ensures robust performance and reproducible results. The concentration of each dNTP is 100 mM.
Brand	Medix Biochemica
Storage	-20°C
Content	1 x 1 ml 100 mM dATP; 1 x 1 ml 100 mM dCTP; 1 x 1 ml 100 mM dGTP; 1 x 1 ml 100 mM dTTP