

Product Datasheet

dUTP (100 mM) MX4503

Product Name	dUTP (100 mM)
Catalog Number	MX4503
Description	MedixMDx 2'-deoxyuridine-5'-triphosphate (dUTP) sodium salt solution is a clear aqueous solution (pH 8.5) with greater than 99% purity. The solution is free of human and bacterial DNA as well as RNases and DNases. This ensures robust performance and reproducible results.
Brand	Medix Biochemica
Storage	-20°C
Content	1 x 1 ml 100 mM dUTP