

Product Datasheet

Anti-SARS-CoV-2 RBD domain Monoclonal Antibody, Azide Free Clone B-B60 857.290

Product Name	Anti-SARS-CoV-2 RBD domain Monoclonal Antibody, Azide Free Clone B-B60
Catalog Number	857.290
Description	Source: Myeloma X63/AG.8653 x Balb/c node cells
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 RBD domain (S1 Spike protein)
Alternative Names	Receptor binding domain Spike protein RBD S1 protein
Brand	Diaclone
Form/Appearance	Liquid
Concentration	1.0 mg/ml (+/- 10 %)
Storage	Stable at +2-8 °C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Buffer	Phosphate-buffered saline. Sterile-filtered through 0.22µm. Carrier and preservative free.
Specificity	RBD domain from Spike protein (aa319-541)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells