

Product Datasheet

Anti-Human CD2 Monoclonal Antibody, Azide Free Clone B-E2 852.500.005

Product Name	Anti-Human CD2 Monoclonal Antibody, Azide Free Clone B-E2
Catalog Number	852.500.005
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	FC, Functional assay, IHC
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Source	Human thymocytes
Specificity	Recognises the Lymphocyte Function associated Antigen-2 (LFA-2), a 50 kDa protein
Biological Activity	Inhibits mixed lymphocyte culture
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells
Immunisation	Human thymocytes