

Product Datasheet

GreenDye 20x 2000

Product Name GreenDye 20x

Catalog Number #2000

Description GreenDye is a highly sensitive intercalating DNA dye with a very low polymerase inhibition. GreenDye features very high fluorescence, consistent results and thermostability for more than 50 cycles. It is compatible with real-time PCR cyclers that support SYBR or FAM channels making it suitable for numerous applications.

For more information, please check :
GreenDye 20x

Brand myPOLS Biotec

Storage -20°C

Content S pack: 1 x 250 µl GreenDye 20x
M pack: 1 x 1000 µl GreenDye 20x