

Product Datasheet

Anti-Human CD2 Monoclonal Antibody, Biotin Conjugated Clone B-E2 852.504.010

Product Name	Anti-Human CD2 Monoclonal Antibody, Biotin Conjugated Clone B-E2
Catalog Number	852.504.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Biotin
Tested Applications	IHC, FC
Immunogen	Human thymocytes
Brand	Diaclone
Form/Appearance	Lyophilised in phosphate-buffered saline with 5% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for six months after reconstitution. For longer storage freeze aliquots.
Note	Reconstitute with 1 ml deionised water Use 10 µl to label 106 cells or 100 µl of whole blood
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Lymphocyte Function associated Antigen-2 (LFA-2), a 50 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells