

Product Datasheet

Anti-Human CD11a Monoclonal Antibody, Azide Free Clone B-B15 852.510.005

Product Name Anti-Human CD11a Monoclonal Antibody, Azide Free Clone B-B15

Catalog Number 852.510.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen Human PBL

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the αL integrin of LFA-1 complex, a 180 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells