

Product Datasheet

RevTaq RT-PCR DNA polymerase #6500

Product Name	RevTaq RT-PCR DNA polymerase
Catalog Number	#6500
Description	<p>RevTaq RT-PCR DNA polymerase is an engineered, extremely thermostable, aptamer-based fast-formulated reverse transcriptase and combined DNA polymerase with a half-life at 95°C of > 40 min. This enzyme allows reverse transcription reactions at high temperatures, minimizing the problems encountered with strong secondary structures in RNA. RevTaq RT-PCR DNA polymerase allows “zero-step” RT-PCRs directly from RNA templates without an isothermal reverse transcription step, as reverse transcription takes place simultaneously with DNA amplification during the cycled PCR elongation step. RevTaq RT-PCR DNA polymerase is the pure, reverse transcription active enzyme ingredient of our Volcano3G® RT-PCR Master Mixes.</p> <p>For more information, please check : RevTaq RT-PCR DNA polymerase, 50x concentrated glycerol stock</p>
Tested Applications	End-Point, Real-Time
Brand	myPOLS Biotec
Storage	-20°C
References	Blatter, N., Bergen, K., Nolte, O., Welte, W., Diederichs, K., Mayer, J., ... & Marx, A. (2013). Structure and Function of an RNA-Reading Thermostable DNA Polymerase. <i>Angewandte Chemie International Edition</i> , 52(45), 11935-11939.
Content	<p>S pack: 1 x 50 µl RevTaq RT-PCR DNA polymerase (50x); 1 x 1.25 ml 5x Volcano buffer</p> <p>L pack: 1 x 250 µ RevTaq RT-PCR DNA polymerase (50x); 2 x 1.25 ml 5x Volcano buffer</p>