

## Product Datasheet

### Phencyclidine (PCP) Antibody 100934

<b>Product Name</b>	Phencyclidine (PCP) Antibody
<b>Catalog Number</b>	100934
<b>Description</b>	Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components
<b>Tested Applications</b>	LF, ELISA
<b>Immunogen</b>	Phencyclidine-BgG
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	≥ 1.5 mg/ml
<b>Storage</b>	2–8 °C
<b>Note</b>	This product was previously available with Catalog code P01-99-65M-P. It has been produced using the same original mouse monoclonal hybridoma cell line.
<b>Isotype</b>	IgG3
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 %
<b>Cross Reactivity</b>	No cross reactivity with Benzoylecgonine, opiates, cannabinoids, amphetamines
<b>Specificity</b>	Antibody primarily recognizes phencyclidine
<b>Matched Pairs</b>	170090 (HRP conjugate as detection)
<b>Shelf Life</b>	36 months
<b>Biological Activity</b>	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN <sub>3</sub> as a preservative