

Product Datasheet

Phencyclidine (PCP) Antibody 100934

Product Name	Phencyclidine (PCP) Antibody
Catalog Number	100934
Description	Protein A purified monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Phencyclidine-BgG
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	≥ 1.5 mg/ml
Storage	+2-8 °C
Note	This product was previously available with Catalog code P01-99-65M-P. It has been produced using the same original mouse monoclonal hybridoma cell line.
Isotype	IgG3
Clonality	Monoclonal
Epitope	Not determined
Purity	≥ 95 %
Buffer	0.05 M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative
Cross Reactivity	No cross reactivity with Benzoyllecgonine, opiates, cannabinoids, amphetamines
Specificity	Antibody primarily recognizes phencyclidine
Matched Pairs	170090 (HRP conjugate as detection)
Shelf Life	36 months
Biological Activity	0.05M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl, 0.05% NaN ₃ as a preservative