

Product Datasheet

Cotinine-4-HRP Antigen 170049

Product Name	Cotinine-4-HRP Antigen
Catalog Number	170049
Description	Synthetic Cotinine, HRP conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-HRP conjugate
Storage	Storage at +2-8 °C. Do not freeze.
Note	This product was previously available with a Catalog code P92-99-02H. It has been produced using the same conjugation protocol.
Buffer	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016 % Gentamicin sulfate, 0.005 % thimerosal
Matched Pairs	PP194IF, PP195IF, PP198IF (capture)
Shelf Life	Unspecified