

Product Datasheet

Cannabinoid (THC)-BSA Antigen 170044

Product Name	Cannabinoid (THC)-BSA Antigen
Catalog Number	170044
Description	Synthetic Cannabinoid, BSA conjugated
Brand	Medix Biochemica
Form/Appearance	Liquid; Antigen-BSA conjugate
Concentration	Lot-specific
Storage	Short term storage (≤ 2 weeks): 2–8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
Note	This product was previously available with a Catalog code P81-99-21A. It has been produced using the same conjugation protocol.
Buffer	Distilled water
Matched Pairs	100892, 100893, 100894, 100898 or P01-99-21M-P (detection)
Shelf Life	60 months