

Product Datasheet

Zika Virus (ZIKV-NS1) Antibody HM1215

Product Name	Zika Virus (ZIKV-NS1) Antibody
Catalog Number	HM1215
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Native Zika Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Dengue, Chikungunya, Yellow Fever, Japanese Encephalitis, Tick-borne Encephalitis, or West Nile Virus
Specificity	Zika Virus NS1
Matched Pairs	HM1216, HM1217, HM1218 or HM1219 (capture)
Preservatives	0.05% Sodium Azide