

Product Datasheet

Digoxin Antibody HM1213

Product Name	Digoxin Antibody
Catalog Number	HM1213
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Digoxin-BSA
Alternative Names	dig
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Specificity	Digoxin
Preservatives	0.05% Sodium Azide