

Product Datasheet

Digoxin Antibody HM1214

| | |
|----------------------------|--|
| Product Name | Digoxin Antibody |
| Catalog Number | HM1214 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | Digoxin-BSA |
| Alternative Names | dig |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | > 95 % by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl |
| Specificity | Digoxin |
| Preservatives | 0.05% Sodium Azide |