

Product Datasheet

Neurofilament H (NfH) Antibody HM1249

| | |
|----------------------------|--|
| Product Name | Neurofilament H (NfH) Antibody |
| Catalog Number | HM1249 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Immunogen | Recombinant NfH protein |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Matched Pairs | HM1254 (detection) |
| Preservatives | 0.05% Sodium Azide |