

# Product Datasheet

## Neurofilament H (NfH) Antibody HM1254

Product Name	Neurofilament H (NfH) Antibody
Catalog Number	HM1254
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant NfH protein
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1249, HM1250, HM1252, HM1253 (capture)
Preservatives	0.05% Sodium Azide