

Product Datasheet

Neurofilament H (NfH) Antibody HM1254

Product Name	Neurofilament H (NfH) Antibody
Catalog Number	HM1254
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant NfH protein
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Matched Pairs	HM1249, HM1250, HM1252, HM1253 (capture)
Preservatives	0.05% Sodium Azide