

Product Datasheet

Neurofilament H (NfH) Antigen, Recombinant >90% LA789

Product Name	Neurofilament H (NfH) Antigen, Recombinant >90%
Catalog Number	LA789
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Epitope	amino acids 93-356
Purity	>85%
Buffer	10 mM Phosphate Buffered Saline with 20% glycerol
Molecular Weight	MW= 62.9 kDa
Matched Pairs	HM1247-HM1254
Expression System	E. coli
Preservatives	0.05% Sodium Azide