

Product Datasheet

Alpha-fetoprotein (AFP) Antibody HM1257

Product Name	Alpha-fetoprotein (AFP) Antibody
Catalog Number	HM1257
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	AFP
Alternative Names	a-fetoprotein, alfa-1-fetoprotein, alpha-fetoglobulin, alpha fetal protein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	AFP
Matched Pairs	HM1258
Preservatives	0.09% Sodium Azide