

Product Datasheet

Prostate Specific Antigen (PSA) Antibody HM1283

Product Name	Prostate Specific Antigen (PSA) Antibody
Catalog Number	HM1283
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA
Immunogen	PSA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Specificity	PSA
Preservatives	0.09% Sodium Azide