

## Product Datasheet

### Procalcitonin (PCT) Antibody, Canine HM1205

<b>Product Name</b>	Procalcitonin (PCT) Antibody, Canine
<b>Catalog Number</b>	HM1205
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM1203, HM1204 or HM1206
<b>Preservatives</b>	0.09% Sodium Azide