

Product Datasheet

Progesterone, Feline Antibody HM1336

Product Name	Progesterone, Feline Antibody
Catalog Number	HM1336
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Progesterone-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	Progesterone
Matched Pairs	LA330 (capture)
Preservatives	0.05% Sodium Azide