

Product Datasheet

Thyroxine (T4), Canine Antibody HM1401

Product Name	Thyroxine (T4), Canine Antibody
Catalog Number	HM1401
Description	Purified Mouse Monoclonal Antibody, in vitro cell culture
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	T4-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Mouse Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	T4
Preservatives	0.1% Sodium Azide