

## Product Datasheet

### Thyroxine (T4), Canine Antibody HM1401

<b>Product Name</b>	Thyroxine (T4), Canine Antibody
<b>Catalog Number</b>	HM1401
<b>Description</b>	Purified Mouse Monoclonal Antibody, in vitro cell culture
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	T4-BSA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Mouse Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 %
<b>Buffer</b>	0.05M Phosphate Buffered Saline, pH 7.2 ± 0.2
<b>Specificity</b>	T4
<b>Preservatives</b>	0.1% Sodium Azide